


ProForm NX™


Modification Guidelines for Professionals in Seating and Positioning

Designed for people requiring greater positioning support and/or cushion modifications. Examples of common modifications:




Wind Swept Deformity

- Cut the Contoured Base to accommodate the pelvic rotation.
- Cut down the existing medial thigh separator.
- Create a new medial thigh separator using Cheat Sheets™ accessory pads.
- Cut Modifiable Thigh Cushion and Modifiable Wedge as needed to conform to the Contoured Base.
- Cinch cover.




Cerebral Vascular Accident (CVA) Above Knee Amputation (AKA) Limited Hip Flexion

- Cut Contoured Base and Modifiable Thigh Cushion up the middle to ischial bar.
- For CVA and AKA, cut the Modifiable Wedge and place only under the affected leg. For Limited Hip Flexion, place Modifiable Wedge only under the unaffected leg.




Symmetrical Posture Requiring Additional Support

- Modify the Contoured Base with Cheat Sheets™ accessory pads to build up the medial thigh separator and/or laterals.



Existing Reddened Skin

- Identify location of bony prominence relative to the base, ideally with a pressure mapping system.
- Cut holes or recessed areas in the Contoured Base to provide additional pressure relief.



Upper Leg Length Discrepancy

- Cut along the Modification Guides on the Contoured Base to accommodate shorter femur.
- Cut the Modifiable Thigh Cushion and Modifiable Wedge to conform to the Contoured Base.
- Cinch cover. (see steps 4-6)



Why select the ProForm NX cushion?

Seating specialists choose the ProForm NX cushion because it can be modified in an additive or subtractive manner. For example, it can provide additional medial or lateral thigh support by building up areas of the Contoured Base with VARILITE Cheat Sheets™ accessory pads. The ProForm NX cushion also can be modified by cutting the Contoured Base, Modifiable Thigh Cushion, and Modifiable Wedge to address a leg length discrepancy.



The ProForm NX cushion is the customizable, out-of-the-box, seating solution.

Modifying Components

1. Modify Contoured Base with a blade or electric knife.
2. Cut Modifiable Wedge with a blade or electric knife.
3. Cut Modifiable Thigh Cushion with a blade or electric knife. Use scored lines on bottom for dimensional reference.



Use care when cutting with a blade or knife.

Cinching the Cover

4. Fold provided hook material in half so that you have hook on both sides.
5. Place hook material on loop at front, bottom edge of cushion on side of modification.
6. Pull firmly on cover and fold it back on itself and attach hook to appropriate loop.

ProForm NX is available with the following options:

- | | |
|------------------------------|--------------------------------|
| Cover: | Mesh or Incontinence |
| Air-Foam Floatation Cushion: | Single-Chamber or Dual-Chamber |

See price list for size availability
HCPCS Code K0736/K0737

